

# Corflo® Enteral Range Update

**BULLETIN**  
May 2008

The NPSA/2007/19 deadline of 31st March 2008 meant that all enteral feeding devices must be:

- Non-compatible with IV syringes and
- Labelled "Enteral"

Since the 1st April, **all** CORFLO enteral feeding devices, and the CE marked pH Indicator strips, supplied by Merck Serono are compatible with NPSA guidelines.

## The CORFLO enteral feeding devices now conform to the following NPSA specifications:

- Not compatible with IV syringes
- Compatible with catheter tip syringes
- Compatible with female luer syringes
- Labelled "Enteral" (and colour coded) ports
- Completely radio-opaque (Fine Bore Feeding Tube range)
- Regular centimetre markings

These features are available across the Merck Serono enteral feeding device range, so you can now order with confidence knowing that you will receive products that will benefit your health care practice by reducing the risk of wrong route errors.



**NPSA Compliant**

**Look out for this symbol! We'll be using it whenever we refer to our range of NPSA compliant products.**

If you require any further information then please contact the local Merck Serono representative whose contact details can be found at [www.MerckGE.co.uk/contact](http://www.MerckGE.co.uk/contact) us or call the

**Merck Serono Support Line FREE on 0800 030 4368**